

Electronic Scale

Leading the Way in Accuracy®

Instruction Manual

Thank you for purchasing a Taylor® Electronic Scale. All of our precision instruments for weight measurement are engineered and designed to meet the highest quality standards ...to assure you years of uncompromising accuracy and consistently dependable, convenient performance. Please read this instruction manual carefully before use. Keep these instructions handy for future reference.

General Operation


1. This scale operates on 3 "AAA" batteries, which are included. Please see "Battery Installation" section of this manual for installation instructions. Also, some models have a static cling label on the lens to prevent scratching. Please remove before use.
2. Your scale has been set at the factory to read weight in pounds. You may select kilogram of measurement by pressing the button on the bottom of the scale.
3. Place scale on flat, hard surface. Carpeted or uneven floors may affect accuracy.
4. Step on the scale (disregarding the reading) and step off. This process insures the scale is set and is now ready for use. Do this whenever the scale is moved. At all other times step straight onto the scale.
5. Step onto the scale platform. The display will count up to your weight. The display will flash twice when the weight is locked. Your weight is now displayed.
6. The scale will turn off automatically. **NOTE:** If you've weighed yourself and a "C" is displayed after the zero, weigh yourself again for a more accurate reading. If the scale displays "Err", this indicates the weight is not stable and the scale cannot lock in the weight. Step off and back onto the scale, standing still while your weight computes.

Low Battery

Your Taylor® scale is equipped with a low battery indicator. It will display "▢" or readings will grow dim or become irregular if the batteries need replacement.

Battery Installation

1. Use 3 AAA batteries (included). Remove any plastic wrap from the batteries before proceeding). Always replace all 3 batteries at the same time; do not combine old and new batteries.
2. Open the battery compartment on the bottom of the scale.
3. Place the new batteries into the battery compartment according to the "+" and "-" polarity symbols at the bottom of the compartment.
4. Replace the battery cover.



Trouble Shooting

In Case of Difficulty...

1. Place your feet in the same position every time you weigh yourself. Stand still during measurement.
2. Weigh yourself without clothing whenever possible. Clothes can vary in weight.
3. Move the scale to a location where it won't be bumped. If your scale is bumped while being activated or while in use, an inaccurate reading may occur.
4. For the most consistent readings, weigh yourself at the same time of day. Allow for weight fluctuations due to food or drink in your system.
5. The maximum capacity of this scale is 440 lbs / 200 kgs. Should "O-LD" appear on the display during a weighing, this indicates that the maximum weighing capacity of the scale has been exceeded. Remove the weight immediately; otherwise, permanent damage may occur.

Reaching Target Weight

Maintaining proper weight is essential to healthy living. Obesity has been linked to a higher risk of heart disease, diabetes and some forms of cancer. The Body Mass Index (BMI) is the most widely accepted measurement of weight according to health professionals.

Use the chart below to determine your own BMI by matching your height in the left hand column with your weight in the center. Anyone, male or female, with a BMI of 25 to 29 is considered over weight and a body mass index of 30 is considered obese. Remember to eat right and exercise.

Body Mass Index

Height (inches)	Body Weight (pounds)						Overweight					Obese					
	91	96	100	105	110	115	119	124	129	134	138	143	148	153	158	162	167
58	91	96	100	105	110	115	119	124	129	134	138	143	148	153	158	162	167
59	94	99	104	109	114	119	124	128	133	138	143	148	153	158	163	168	173
60	97	102	107	112	118	123	128	133	138	143	148	153	158	163	168	174	179
61	100	106	111	116	122	127	132	137	143	148	153	158	164	169	174	180	185
62	104	109	115	120	126	131	136	142	147	153	158	164	169	175	180	186	191
63	107	113	118	124	130	135	141	146	152	158	163	169	175	180	186	191	197
64	110	116	122	128	134	140	145	151	157	163	169	174	180	186	192	197	204
65	114	120	126	132	138	144	150	156	162	168	174	180	186	192	198	204	210
66	118	124	130	136	142	148	155	161	167	173	179	186	192	198	204	210	216
67	121	127	134	140	146	153	159	166	172	178	185	191	198	204	211	217	223
68	125	131	138	144	151	158	164	171	177	184	190	197	203	210	216	223	230
69	128	135	142	149	155	162	169	176	182	189	196	203	209	216	223	230	236
70	132	139	146	153	160	167	174	181	188	195	202	209	216	222	229	236	243
71	136	143	150	157	165	172	179	186	193	200	208	215	222	229	236	243	250
72	140	147	154	162	169	177	184	191	199	206	213	221	228	235	242	250	258
73	144	151	159	166	174	182	189	197	204	212	219	227	235	242	250	257	265
74	148	155	163	171	179	186	194	202	210	218	225	233	241	249	256	264	272
75	152	160	168	176	184	192	200	208	216	224	232	240	248	256	264	272	279
76	156	164	172	180	189	197	205	213	221	230	238	246	254	263	271	279	287
	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35

Body Mass Index

Precautions

1. Clean exposed parts with a soft, slightly damp cloth. To remove stains, use a mild soap. Never use detergents, excess water, treated cloths, harsh cleaning agents, or sprays. Do not attempt to lubricate or open scale casing as this will void your warranty.
2. Although your scale is designed to be maintenance-free, it contains sensitive electronic parts. Avoid rough treatment. Do not drop, kick, or jump on it. Treat it with care to ensure the best performance.
3. This scale is designed to read weights between 25 and 440 lbs / 12 and 200 kgs. To prevent damage, do not place more than 440 lbs / 200 kg on it.
4. Do not store scale where you store cleaning chemicals. The vapors from some household products may affect the electronic components of your scale. Do not store scale on its side.
5. This scale is an extremely sensitive weighing device. In order to prevent run down of the batteries, do not store anything on the scale.
6. Do not dispose of batteries in fire. Batteries may explode or leak. Do not mix Alkaline, carbon-zinc (standard) or Nickel-Cadmium (rechargeable) batteries. Remove the batteries if the scale will not be used for a long period of time.

Five (5) Year Limited Warranty

This scale is warranted against defects in materials or workmanship (excluding batteries) for five (5) years of the original purchaser from date of retail purchase. It does not cover damages or wear resulting from accident, misuse, abuse, commercial use, or unauthorized adjustment and/or repair. Do not return to retailer. Should this scale require service (or replacement at our option) while under warranty, please pack the item in the original packaging and return it prepaid, along with store receipt showing date of purchase and a note explaining reason for return to:

TAYLOR Precision Products
2220 Entrada Del Sol, Suite A
Las Cruces, New Mexico 88001 USA

There are no express warranties except as listed above. This warranty gives you specific legal rights, and you may have other rights which vary from state to state. For additional product information, or warranty information outside the USA, please contact us through www.taylorusa.com.

©2011 Taylor Precision Products and its affiliated companies, all rights reserved. Taylor® and Leading the Way in Accuracy® are registered trademarks of Taylor Precision Products and its affiliated companies. All rights reserved.

Made to our exact specifications in China.
Not legal for trade.
www.taylorusa.com